

PROSTHETIC SOCKS

KCN | KCM | KCW



The ALPS® prosthetic sock is lightweight and soft, providing cushioning and comfort for a superior fit.

PRODUCT SPECIFICATIONS

NARROW, MEDIUM, & WIDE

1, 3, & 5 PLY


EXTRA SHORT, SHORT, REGULAR, & LONG

THREE STYLES

THREE SIZES

THREE THICKNESSES

FOUR LENGTHS



DISTAL HOLE (H)




REINFORCED DISTAL HOLE (HR)



NO HOLE (N)

HIGHLIGHTS

-  Assists with volume fluctuation
-  Technical fibers allow a cool, comfortable wear
-  Keeps its shape over time; minimal stretch and will not mat down
-  **NOW AVAILABLE IN BUNDLES!**
A complete solution to patients to manage the volume fluctuation of their residual limb. Most popular thicknesses are available on hand.



MADE WITH PREMIUM TECHNICAL FIBERS

CREATES A MOISTURE MANAGEMENT SYSTEM TO MOVE PERSPIRATION AWAY FROM THE BODY AND THROUGH THE FABRIC, WHERE IT CAN EVAPORATE QUICKLY.

WITH BUNDLES, YOUR PATIENT CAN:


- Select the socks needed to manage their variation
- Have prosthetic socks always available when others are washing



making lives better™


FITTING DIAGRAM

It's recommended the sock to be approximately 5 cm longer than the socket



HOW TO SIZE

Measurements taken 5cm from the distal end



SIZING CHART

FIND THE PERFECT FIT.

RECOMMENDED ALPS LINER SIZE	LEVEL OF AMPUTATION	PROXIMAL CIRCUMFERENCE	DISTAL CIRCUMFERENCE AT (5 CM)	*LENGHT OF THE SOCKET 13-18 CM	*LENGHT OF THE SOCKET 18-23 CM	*LENGHT OF THE SOCKET 23-38 CM	*LENGHT OF THE SOCKET 38-48 CM
20	BK	min. 20-max. 40 cm	20-23 cm	KCN-8-X-Y	KCN-10-X-Y	KCN-14-X-Y	KCN-18-X-Y
24, 26, 28	BK	min. 30-max. 60 cm	23-31 cm	KCM-8-X-Y	KCM-10-X-Y	KCM-14-X-Y	KCM-18-X-Y
32, 38	AK	min. 40-max. 80 cm	31-36 cm	KCW-8-X-Y	KCW-10-X-Y	KCW-14-X-Y	KCW-18-X-Y

X indicates the thickness in Ply (1Ply= 1,5mm; 3Ply= 2mm; 5Ply=2,5mm)

Y indicates the distal end (N-without hole; H with hole; HR with reinforced hole)

EXAMPLES:

KCN-08-3-H Prosthetic sock, Narrow, Extra Short, 3 Ply, with hole

KCM-10-1-N Prosthetic sock, Medium, Short, 1 Ply, without hole

KCW-18-5-HR Prosthetic Sock, Wide, Long, 5 Ply, with reinforced hole



Tel +420 377 223 127 | Fax +420 377 223 010

www.easyliner.eu | info@easyliner.eu

Božkovské náměstí 17/21, 32600 Plzeň – Czech Republic

