

EVERYDAY

Freedom DynAdapt™

Fit, Finish and Roll

- Slim profile, carbon fiber foot system with a slender, anatomic design
- Full length attached heel with no bolts or bushings
- Combination of smooth rollover function, improved stability, and dynamic response makes DynAdapt the best all-around foot for all K3 users
- EnduraCore® Technology delivers up to 3 times the fatigue life of a standard laminate
- Split keel and heel design delivers enhanced stability
- Lightweight construction – only 635 grams (incl. foot shell)
- Rated for 166 kg
- 36-month warranty



Freedom DynAdapt (F10)



Split Toe




Sandal Toe



EnduraCore®
Technology

Freedom DynAdapt™



Note: Build height measurements represent a fully assembled unit.

Product ID: F10

Sizes: 22-31 cm

Sandal toe available for sizes 22-28 cm*

User Weight Rating: 166 kg

Heel Height: 10 mm

Connector: Male pyramid

Warranty: Foot module 36 months,
Cosmetic foot shell 6 months

*When ordering the Sandal Toe option, please contact Customer Service for lead time.

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
 - 22-25 cm (KIT-00-11373)
 - 26-28 cm (KIT-00-11374)
 - 29-31 cm (KIT-00-11375)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

| Size (cm) | Height (mm) | Weight (g) |
|-----------|-------------|------------|
| 22-25 | 158 | 535 |
| 26-28 | 170 | 635 |
| 29-31 | 182 | 750 |

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

| Impact Level | | | Weight (kg) |
|--------------|----------|------|-------------|
| Low | Moderate | High | |
| 1 | 1 | 2 | 44-52 |
| 1 | 2 | 3 | 53-59 |
| 2 | 3 | 4 | 60-68 |
| 3 | 4 | 5 | 69-77 |
| 4 | 5 | 6 | 78-88 |
| 5 | 6 | 7 | 89-100 |
| 6 | 7 | 8 | 101-116 |
| 7 | 8 | 9 | 117-130 |
| 8 | 9 | | 131-147 |
| 9 | | | 148-166 |

Complete Prosthetic Foot (incl. foot shell) Part Number

F10-00--A---

Toe: 0-regular, 1-sandal
 Category: 1-9
 Size: 22-31
 Version: M-mid profile, S-sandal toe
 Side: L-left, R-right
 Generation: 3-regular toe, 2-sandal toe
 Skin Tone: L-light, M-medium, D-dark

Foot Module Part Number

F10-00-0A--

Category: 1-9
 Size: 22-31
 Toe Configuration: RU-regular toe, universal, SR-sandal toe, right, SL-sandal toe, left

Foot Shell Part Number

FTC----4--

Generation: 2-3
 Profile: M-mid, H-high
 Cap: 0-yes, 1-no
 Size: 22-31
 Skin Tone: L-light, M-medium, D-dark
 Toe: R-regular, S-sandal
 Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.