


Realistic Prosthetics Ltd.

Dimension Chart for RealSkin PVC Myo / Mech. Cosmetic Gloves

GLOVE	Otto Bock Size *	Hand	Across Knuckles	A	B	C	D
				Palm Circumference	Wrist Circumference	Middle Finger	Length to Thumb Tip
RPL 401	8 ¼	Left	3 ½" 89 mm	219 mm	176 mm	76 mm	580 mm
RPL 402	8 ¼	Right	3 ½" 89 mm	219 mm	176 mm	76 mm	580 mm
RPL 403	7 ¾	Left	3 ¼" 82.5 mm	205 mm	168 mm	71 mm	570 mm
RPL 404	7 ¾	Right	3 ¼" 82.5 mm	205 mm	168 mm	71 mm	570 mm
RPL 405	7 ¼	Left	3 " 76 mm	194 mm	168 mm	71 mm	540 mm
RPL 406	7 ¼	Right	3 " 76 mm	194 mm	168 mm	71 mm	540 mm
RPL 407	6 ¾	Left	2 ¾" 70 mm	174 mm	150 mm	67 mm	440 mm
RPL 408	6 ¾	Right	2 ¾" 70 mm	174 mm	150 mm	67 mm	440 mm

* Otto Bock is the Trade Mark of Otto Bock Healthcare Plc.



Our cosmetic gloves are compatible with UK, European and other manufacturers hand systems
All our gloves are CE marked

Realistic Prosthetics Ltd.

Unit 19 Platts Eyot, Lower Sunbury Road, Hampton, Middlesex TW12 2HF U.K.

Phone: 020 8783 0099

Fax: 020 8941 4077

International +44 (0)

Web: www.realisticprosthetics.com

E-mail: info@realisticprosthetics.com